



ProFem[®]

PELVIC FLOOR CARE - BY SAYCO



Cube

With Drains





Cube

Without Drains





The Cube Pessary, available with or without drains, may be used to relieve the symptoms of a cystocele, or rectocele as well as a second or third-degree uterine prolapse. Your health care professional will teach you self care.

INSERTION 1-3



REMOVAL 4-6



This Cube should be removed, cleaned and re-inserted daily, to maintain hygienic conditions.

All white silicone pessaries in this range are for single patient use.

mm	Code No.	mm	Code No.
25	CU25D #0	25	CU25 #0
29	CU29D #1	29	CU29 #1
33	CU33D #2	33	CU33 #2
37	CU37D #3	37	CU37 #3
41	CU41D #4	41	CU41 #4
45	CU45D #5	45	CU45 #5
50	CU50D #6	50	CU50 #6
56	CU56D #7	56	CU56 #7
63	CU63D #8	63	CU63 #8
70	CU70D #9	70	CU70 #9
75	CU75D #10	75	CU75 #10



Sayco Pty Ltd Unit 5, 260 Wickham Road, Moorabbin, Victoria 3189, Australia
Toll Free: 1800 333 456 Email: customerservice@profem.com.au

www.profem.com.au