

# ProFem<sup>®</sup>

PELVIC FLOOR CARE - BY SAYCO



## Ring

Without Support



## Ring

With Support




The Ring pessary, available with or without a supporting membrane, is primarily used for a mild first-degree uterine prolapse or a cystocele.

The Profem Ring pessary has two notches to allow for the pessary to be easily folded.

mm	Inch	Code No.	mm	Inch	Code No.
43	1.75	R1.75 #0	43	1.75	R1.75S #0
51	2.00	R2.00 #1	51	2.00	R2.00S #1
56	2.25	R2.25 #2	56	2.25	R2.25S #2
62	2.50	R2.50 #3	62	2.50	R2.50S #3
69	2.75	R2.75 #4	69	2.75	R2.75S #4
74	3.00	R3.00 #5	74	3.00	R3.00S #5
81	3.25	R3.25 #6	81	3.25	R3.25S #6
87	3.50	R3.50 #7	87	3.50	R3.50S #7
95	3.75	R3.75 #8	95	3.75	R3.75S #8
100	4.00	R4.00 #9	100	4.00	R4.00S #9

All white silicone pessaries in this range are for single patient use.



Sayco Pty Ltd Unit 5, 260 Wickham Road, Moorabbin, Victoria 3189, Australia  
Toll Free: 1800 333 456 Email: customerservice@profem.com.au

www.profem.com.au